


The TAD's are inserted directly through the gingiva using either topical or local anesthetic. The procedure is absolutely pain free. No incisions are necessary and therefore no sutures. Occasionally, patients will have a feeling of pressure.



After treatment the titanium TAD's can easily be removed without anesthetic.

YOUR NEXT APPOINTMENT DATE

	Date	Time
MON		
TUE		
WED		
THUR		
FRI		

The TAD's will be securely anchored to the bone at the end of the procedure. To ensure success and prevent loosening, the patient must pay close attention to two things:

1. Excellent Oral Hygiene

Please clean the area around the TAD thoroughly using a toothbrush and toothpaste (just like you would a tooth) in order to prevent the accumulation of plaque and bacteria.

2. Do not disturb the TAD

Please do not move the TAD by playing with it with your tongue or fingers. The TAD is more likely to come loose if it is disturbed.

Your specialist for TAD's:

Tannlege Matthias Möller, Dr. med. dent.
Spesialist i kjeveortopedi

Tannhelsetjenestens kompetansesenter Sør
Arendal, Norway

More information:

Distributed by:

Manufactured by:



Moltkestraße 41
78532 Tuttlingen, Germany
Phone +49 7461 96637-0
Fax +49 7461 96637-29
E-Mail info@psm.ms
Internet www.psm.ms

©FRICOM.de 12.10

TAD's (Temporary Anchorage Devices)
MICRO IMPLANTS



- Reduce the need for extractions
- Reduce the need for headgear
- Create beautiful, symmetrical smiles
- Work well for children and adults

BENEFIT / LOMAS


*The most innovative anchorage method
in orthodontics!*

The most important principle in orthodontics is anchorage

When teeth or groups of teeth have to be moved, orthodontists need solid anchorage.

Creating reliable anchorage is one of the most important things orthodontists do.

Orthodontic devices like headgears or face masks can be used as anchorage sources. Patients, often times, find them annoying, and therefore do not wear them sufficiently, preventing successful achievement of treatment goals. In these situations often times the alternative is the extraction of healthy permanent teeth. TAD's now provide a new and better alternative for both of these issues.



Oftentimes, headgears and face masks are no longer needed thanks to TAD's!

Advantages:

- Invisible when inserted
- Reduces or eliminates the need for headgear
- Reduced treatment time
- Easy oral hygiene
- Reduced need for extractions
- Safe and easy insertion
- Successful long term history
- Made in Germany

Innovative and safe with the Benefit System

Benefit TAD's are ideal to solve anchorage requirements in an easy and **invisible** way. Unlike dental implants, where you must wait several months or more before they can be used, titanium TAD's can be used immediately after insertion. This is because they do not osseointegrate and this also makes them easy to remove.

The Benefit System consists of titanium TAD's and various abutments which can be attached to them. For extra stability, two TAD's can be used which are connected using the BENEplate.



BENEFIT Abutments



1:1

BENEFIT Micro-Implants



BENEPLATE



1:1

LOMAS Micro-Implants

BENEFIT in various situations

Gaining space without tooth extraction: The BENESlider



before



after

"Invisible" space closure: The MESIALslider



before



after 10 months

